

Equine wound care: The need for a practical solution

Conventional Treatments: Bandages, Ointments and Creams

Re-application/re-bandaging

- Constant re-application and re-bandaging
- Owner/patient compliance with treatment regimen

Difficult to bandage areas

- Wraps and bandages are ill suited for areas such as the coronary band, hocks and other joints



“The ideal wound dressing would function like skin; to protect the area from trauma, thermoregulate, allow water and electrolytes to pass through, and allow oxygen and carbon dioxide exchange”

Paul Wan, DVM, DACVS

Equine Veterinary Specialists, Norco, CA




SutureSeal

Seal. Protect. Heal.

For equine wounds and post-op healing

1 application = Up to 14 days of wound protection

 **SutureSeal**
(Hydrogel Wound sealant)



SutureSeal is a hydrogel polymer that protects incisions and wounds for up to 14 days

- Safety – Non-toxic and Biodegradable
- Protection – Flexible, Transparent and Durable (lasts up to 2 weeks)
- Healing – Permits airflow and Locks out contaminants

Before



After 14 days of

 **SutureSeal**

